

rno-miR-29b-3p inhibitor

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	rno-miR-29b-3p inhibitor is a fully methoxy-modified complementary single strand of mature miRNA. This miRNA inhibitor specifically binds mature miRNA, blocking its complementary pairing with target genes, thereby inhibiting miRNA function. It is used for miRNA loss-of-function research.
Targets(IC50)	MicroRNA

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481