

PG-911

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Keep away from moisture

Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

| | |
|---------------|---|
| Description | PG-911 is a full agonist of MC5R with an EC50 of 0.031 nM for hMC5R. It also acts as a complete antagonist of hMC3R and hMC4R, with Kb values of 0.38 nM and 0.17 nM, respectively. PG-911 is applicable in neurological disorder research. |
| Targets(IC50) | Melanocortin Receptor |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481