

CRBN ligand-893

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from direct sunlight
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	CRBN ligand-893 (Compound 4) is a CRBN ligand with a dissociation constant (Kd) of 0.2 μ M. It connects to a target protein ligand via a linker to form a PROTAC molecule and is also applicable for research into molecular glue.
Targets(IC50)	Ligands for E3 Ligase, Molecular Glues, PROTACs

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481