

Glyceraldehyde-3-phosphate Dehydrogenase, Rabbit

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Glyceraldehyde-3-phosphate Dehydrogenase, Rabbit (EC 1.2.1.12), facilitates the conversion of glyceraldehyde-3-phosphate to 1,3-bisphosphoglycerate, which is an integral step in the glycolysis pathway.
Targets(IC50)	Endogenous Metabolite

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481