

Aminopeptidase M, Porcine

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	AminopeptidaseM, Porcine (EC 3.4.11.2), is a metalloprotease that can hydrolyze almost all N-terminal amino acids of unsubstituted oligopeptides. It does not cleave X-Pro bonds or amino acids with blocked N-terminals. AminopeptidaseM is utilized for peptide sequence analysis.
Targets(IC50)	Aminopeptidase

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481