

CDK12/13 ligand 3

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	CDK12/13 ligand 3 is a CDK12/13 PROTAC ligand. It can bind with an E3 ligase ligand and a linker for the synthesis of the CDK12/13 PROTAC degrader DN1679. CDK12/13 ligand 3 is applicable in cancer research.
Targets(IC50)	CDK,Ligands for Target Protein for PROTAC

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481