

4,4'-(Diazoamino)dibzenesulfonic acid disodium

Chemical Properties

CAS No. :	56120-28-6
Formula:	C ₁₂ H ₉ N ₃ Na ₂ O ₆ S ₂
Molecular Weight:	401.32
Storage:	Keep away from direct sunlight Store at -20°C <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	4,4'-(Diazoamino) dibzenesulfonic acid disodium is a compound found in certain colorants.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.4918 mL	12.4589 mL	24.9178 mL
5 mM	0.4984 mL	2.4918 mL	4.9836 mL
10 mM	0.2492 mL	1.2459 mL	2.4918 mL
50 mM	0.0498 mL	0.2492 mL	0.4984 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481