

Mal-PEG-NHS ester (MW 2000)

Chemical Properties

CAS No. : 174459-59-7

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Mal-PEG-NHS ester (MW 2000) is a PEG linker utilized in conjugating lipoteichoic acid (LTA) with CpG DNA fragments to synthesize CpG-TA heterodimers. It is applicable in immune stimulation research.
Targets(IC50)	ADC Linker

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481