

## Celangulin III

## Chemical Properties

CAS No. :	144379-42-0
Formula:	C <sub>34</sub> H <sub>44</sub> O <sub>14</sub>
Molecular Weight:	676.71
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

## Biological Description

Description	Celangulin III is a naturally occurring compound found in <i>Tripterygium wilfordii</i> and <i>Begonia fimbristipula</i> .
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.4777 mL	7.3887 mL	14.7774 mL
5 mM	0.2955 mL	1.4777 mL	2.9555 mL
10 mM	0.1478 mL	0.7389 mL	1.4777 mL
50 mM	0.0296 mL	0.1478 mL	0.2955 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481