

(S)-3-Hydroxytetradecanoyl-CoA

Chemical Properties

CAS No. :	950902-30-4
Formula:	C35H62N7O18P3S
Molecular Weight:	993.89
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	(S)-3-Hydroxytetradecanoyl-CoA is a nucleotide metabolite.
Targets(IC50)	Nucleoside Antimetabolite/Analog

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.0061 mL	5.0307 mL	10.0615 mL
5 mM	0.2012 mL	1.0061 mL	2.0123 mL
10 mM	0.1006 mL	0.5031 mL	1.0061 mL
50 mM	0.0201 mL	0.1006 mL	0.2012 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481