

BCR-ABL degrader 1

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	BCR-ABLdegrader 1 (Compound 23) is a HyT-based BCR-ABL degrader with a DC50 of 19.9 μ M. It is useful for research in chronic myeloid leukemia.
Targets(IC50)	Bcr-Abl

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481