

2-Methylbutanal-d9

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	2-Methylbutanal-d9 (2-Methylbutyraldehyde-d9) is the deuterium-labeled form of 2-Methylbutanal. 2-Methylbutanal (2-Methylbutyraldehyde) serves as a distinctive volatile biomarker, useful for detecting urinary tract pathogens such as Proteus.
Targets(IC50)	Endogenous Metabolite

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481