

## Beclomethasone 21-propionate-d3-1

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Beclomethasone 21-propionate-d3-1 is the deuterium-labeled form of Beclomethasone 21-propionate. This compound is used as a biochemical reagent in glycobiology research, which examines the structure, synthesis, biology, and evolution of sugars. It involves carbohydrate chemistry, glycan formation and degradation enzymology, protein-glycan interactions, and the role of glycans in biological systems, closely connecting with fundamental research, biomedicine, and biotechnology.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481