

## N-Formimino-L-aspartate

## Chemical Properties

CAS No. :	2374-41-6
Formula:	C <sub>5</sub> H <sub>8</sub> N <sub>2</sub> O <sub>4</sub>
Molecular Weight:	160.13
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

## Biological Description

Description	N-Formimino-L-aspartate (N-Formimidoyl-L-aspartate) is a nucleoside metabolite.
Targets(IC50)	Nucleoside Antimetabolite/Analog

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	6.2449 mL	31.2246 mL	62.4493 mL
5 mM	1.249 mL	6.2449 mL	12.4899 mL
10 mM	0.6245 mL	3.1225 mL	6.2449 mL
50 mM	0.1249 mL	0.6245 mL	1.249 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481