

Ro 25-1553 TFA

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Keep away from moisture

Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description

Ro 25-1553 TFA is a 31-amino acid analog of vasoactive intestinal peptide (VIP) that acts as an agonist for the VIP2 receptor (VPAC2 receptor). It exhibits bronchodilatory effects on tracheal smooth muscle contractions induced by neural stimuli or Carbachol in guinea pig models.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481