

## $\alpha$ -Lytic Protease M190A Mutant

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	$\alpha$ -Lytic Protease serves as an alternative specific protease for proteomics applications, with its wild-type (WT) variant cleaving after T, A, S, and V residues. The $\alpha$ -Lytic Protease M190A Mutant, on the other hand, exhibits distinct cleavage specificity, targeting M, F, and L residues.
Targets(IC50)	Endogenous Metabolite

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481