

## Docosanoyl-CoA

## Chemical Properties

CAS No. :	24330-89-0
Formula:	C43H78N7O17P3S
Molecular Weight:	1090.11
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

## Biological Description

Description	Docosanoyl-CoA is a nucleoside metabolite.
Targets(IC50)	Nucleoside Antimetabolite/Analog

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.9173 mL	4.5867 mL	9.1734 mL
5 mM	0.1835 mL	0.9173 mL	1.8347 mL
10 mM	0.0917 mL	0.4587 mL	0.9173 mL
50 mM	0.0183 mL	0.0917 mL	0.1835 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481