

YLRPEGDW

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Keep away from moisture

Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	YLRPEGDW is an orally active peptide that promotes sleep. It can upregulate genes related to autophagy and enhance autophagic activity. Additionally, YLRPEGDW has the ability to counteract sleep reduction caused by Caffeine.
Targets(IC50)	Autophagy

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481