

Cas12a mRNA

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description

Cas12a mRNA encodes the Cas12a protein, which functions as an endonuclease.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481