

## D-4-Hydroxyphenylglycine-d5-1

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	D-4-Hydroxyphenylglycine-d5-1 (D-(-)-4-Hydroxyphenylglycine-d5-1) is a deuterium-labeled form of D-4-Hydroxyphenylglycine. This compound, D-4-Hydroxyphenylglycine (D-(-)-4-Hydroxyphenylglycine), serves as a crucial raw material in the production of semi-synthetic $\beta$ -lactam antibiotics, such as amoxicillin and cefadroxil.
Targets(IC50)	Antibiotic

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481