

Dopamine β -monooxygenase, Human

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Dopamine β -monooxygenase, Human (EC 1.14.17.1), is a copper-containing monooxygenase responsible for converting dopamine into norepinephrine within noradrenergic neurons, adrenergic neurons, and adrenal chromaffin cells.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481