

## 2-Amino-3-methoxybenzoic acid-d3

## Chemical Properties

CAS No. :	2377028-36-7
Formula:	C <sub>8</sub> H <sub>9</sub> NO <sub>3</sub>
Molecular Weight:	170.18
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

## Biological Description

Description	2-Amino-3-methoxybenzoic acid-d3 is a deuterium-labeled version of 2-Amino-3-methoxybenzoic acid. This biochemical reagent can be utilized as a biomaterial or organic compound in research related to life sciences.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.8761 mL	29.3807 mL	58.7613 mL
5 mM	1.1752 mL	5.8761 mL	11.7523 mL
10 mM	0.5876 mL	2.9381 mL	5.8761 mL
50 mM	0.1175 mL	0.5876 mL	1.1752 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481