

## Pulvinamide

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Pulvinamide, a lichen alkaloid pigment identified in <i>P. crocata</i> , exhibits antibacterial and antifungal properties. It is effective against <i>B. subtilis</i> and <i>C. albicans</i> , with minimum inhibitory concentrations (MICs) of 50 µg/ml and 200 µg/ml, respectively.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481