

IL-22 Protein, Mouse, Recombinant (HEK293, His)

General Information

Protein Construction: Leu34-Val179

Species: Mouse

Expression Host: HEK293 Cells

Accession: Q9JJY9

Molecular Weight: 17.73 kDa (Predicted); 30-45 kDa (Reducing condition, due to glycosylation)

QC Testing

Biological Activity: Immobilized Mouse IL-22, His Tag at 2 µg/ml (100 µl/well) on the plate. Dose response curve for Mouse IL-22R alpha 1, hFc Tag with the EC50 of 6.7 ng/ml determined by ELISA (QC Test).

Purity: > 95% as determined by Bis-Tris PAGE; > 95% as determined by HPLC

Endotoxin: < 1 EU/µg of the protein as determined by the LAL method.

Formulation: Lyophilized from 0.22 µm filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.

Preparation and Storage

Reconstitution:

Reconstitute the lyophilized protein in distilled water. The product concentration should not be less than 100 µg/ml. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.

Stability & Storage:

It is recommended to store recombinant proteins at -20°C to -80°C for future use. Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Protein Background

The cytokine interleukin-22 (IL-22) is a critical regulator of epithelial homeostasis. It has been implicated in multiple aspects of epithelial barrier function, including regulation of epithelial cell growth and permeability, production of mucus and antimicrobial proteins (AMPs), and complement production.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481