

## Kallikrein 2/KLK2 Protein, Human, Recombinant

### General Information

Synonyms:	Kallikrein2;KLK-2;Kallikrein-2;hGK-1;KLK2;Glandular Kallikrein-1;KLK 2;Tissue kallikrein2;Glandular kallikrein-1 (hGK-1);Tissue Kallikrein-2;Tissue kallikrein-2
Protein Construction:	Ile25-Pro261
Species:	Human
Expression Host:	HEK293 Cells
Accession:	P20151-1
Molecular Weight:	26.15 kDa (Predicted); 28-35 kDa (Reducing condition, due to glycosylation)

### QC Testing

Biological Activity:	Immobilized Human Kallikrein 2, No Tag at 0.5 µg/ml (100 µl/well) on the plate. Dose response curve for Anti-Kallikrein 2 Antibody, hFc Tag with the EC50 of 3.3 ng/ml determined by ELISA (QC Test).
Purity:	> 95% as determined by Bis-Tris PAGE
Endotoxin:	< 0.1 EU/µg of the protein as determined by the LAL method.
Formulation:	Lyophilized from 0.22 µm filtered solution in 20 mM Tris, 150 mM NaCl, 8% trehalose, 0.05% Brij-35, 5 mM Benzamidine (pH 8.0).

### Preparation and Storage

#### Reconstitution:

Reconstitute the lyophilized protein in distilled water. The product concentration should not be less than 100 µg/ml. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.

#### Stability & Storage:

It is recommended to store recombinant proteins at -20°C to -80°C for future use. Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

#### Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

### Protein Background

Kallikrein-related peptidase 2 (KLK2) is a serine protease exhibiting antiangiogenic properties through proteolytic activity. KLK2 is overexpressed in prostate cancer and plays a pivotal role in cancer progression, establishing it as a potential therapeutic target.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481