

## HLA-A\*02:01&B2M&KRAS WT (KLVVVGAGGV) Monomer Protein, Human, MHC (His & Avi)

### General Information

Synonyms:	K-Ras 2;C-K-RAS;K-RAS4A;RASK2;GTPase Kras;NS;KRAS;CFC2;K-RAS2B;RALD;KI-RAS;NS3;KRAS2;K-RAS4B;K-RAS2A;KRAS1;MHC
Protein Construction:	Gly25-Thr305 (HLA-A*02:01), Ile21-Met119 (B2M) and KLVVVGAGGV peptide
Species:	Human
Expression Host:	E. coli
Accession:	A0A140T913(HLA-A*02:01)&P61769(B2M)&KLVVVGAGGV
Molecular Weight:	35.6 kDa (HLA-A*02:01) and 11.9 kDa (B2M)

### QC Testing

Biological Activity:	Immobilized Human HLA-A*02:01&B2M&KRAS WT (KLVVVGAGGV) Monomer, His Tag at 0.5 µg/ml (100 µl/well) on the plate. Dose response curve for Anti-B2M Antibody, mFc Tag with the EC50 of 18.7 ng/ml determined by ELISA.
Purity:	> 95% as determined by Bis-Tris PAGE; > 95% as determined by HPLC
Endotoxin:	< 1 EU/µg of the protein as determined by the LAL method.
Formulation:	Supplied as 0.22 µm filtered solution in 20 mM Tris, 200 mM NaCl (pH 8.0).

### Preparation and Storage

#### Stability & Storage:

It is recommended to store recombinant proteins at -20°C to -80°C for future use. Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

#### Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

### Protein Background

Kirsten rat sarcoma 2 viral oncogene homolog (KRAS) is the most commonly mutated oncogene in human cancer. The developments of many cancers depend on sustained expression and signaling of KRAS, which makes KRAS a high-priority therapeutic target.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

**This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use**

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481