

HLA-A*02:01&B2M&RAB7B (KLIIVGAIGV) Monomer Protein, Human, MHC (His & Avi), Biotinylated

General Information

Synonyms: RAB7B;MHC

Protein Construction: Gly25-Thr305 (HLA-A*02:01), Ile21-Met119 (B2M) and KLIIVGAIGV peptide

Species: Human

Expression Host: E. coli

Accession: A0A140T913(HLA-A*02:01)&P61769(B2M)&KLIIVGAIGV

Molecular Weight: 35.6 kDa (HLA-A*02:01) and 11.9 kDa (B2M)

QC Testing

Biological Activity: Activity has not been tested. It is theoretically active, but we cannot guarantee it.

Purity: > 95% as determined by Bis-Tris PAGE; > 95% as determined by HPLC

Endotoxin: < 1 EU/μg of the protein as determined by the LAL method.

Formulation: Supplied as 0.22 μm filtered solution in PBS (pH 7.4).

Preparation and Storage

Stability & Storage:

It is recommended to store recombinant proteins at -20°C to -80°C for future use. Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Protein Background

Rab proteins are key-regulators of intracellular membrane trafficking. Rab7b is a recently identified Rab protein that may downregulate TLR4 and TLR9-mediated inflammatory responses. Rab7b, believed to have similar function as Rab7, controls however vesicular trafficking from endosomes to the TGN. It is localized to late endosomes/lysosomes as well as the TGN. Rab7b interferes with enzymes delivery to lysosomes and with the retrograde Shiga toxin transport to the Golgi.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481