

HLA-A*02:01&B2M&P53 R175H (HMTEVVRHC) Tetramer Protein, Human, MHC (His & Avi), APC-

General Information

Synonyms: P53;MHC;TRP53;BCC7;HLA-A;TP53;Antigen NY-CO-13;FLJ92943;LFS1

Protein Construction: Gly25-Thr305 (HLA-A*02:01), Ile21-Met119 (B2M) and HMTEVVRHC peptide

Species: Human

Expression Host: E. coli

Accession: A0A140T913(HLA-A*02:01)&P61769(B2M)&HMTEVVRHC

QC Testing

Biological Activity: Immobilized Anti-HLA-A*02:01&B2M&P53 R175H Antibody, hFc Tag at 2 µg/ml (100 µl/well) on the plate. Dose response curve for APC-Labeled Human HLA-A*02:01&B2M&P53 R175H (HMTEVVRHC) Tetramer, His Tag with the EC50 of 21.1 ng/ml determined by ELISA (QC Test).

Formulation: Supplied as 0.22 µm filtered solution in 20 mM Tris, 200 mM NaCl, 0.2% BSA (pH 8.0).

Preparation and Storage**Stability & Storage:**

It is recommended to store recombinant proteins at -20°C to -80°C for future use. Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Protein Background

p53 is a tumor suppressor protein. Under stressful conditions, p53 tightly regulates cell growth by promoting apoptosis and DNA repair. When p53 becomes mutated, it loses its function, resulting in abnormal cell proliferation and tumor progression. Depending on the p53 mutation, it has been shown to form aggregates leading to negative gain of function of the protein. p53 mutant associated aggregation has been observed in several cancer tissues and has been shown to promote tumor growth.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481