

HLA-A*02:01&B2M&PLK1 (AAPGAPAAA) Monomer Protein, Human, MHC (His & Avi)

General Information

Synonyms:	PLK-1;PLK;STPK13
Protein Construction:	Gly25-Thr305 (HLA-A*02:01), Ile21-Met119 (B2M) and AAPGAPAAA peptide
Species:	Human
Expression Host:	E. coli
Accession:	A0A140T913(HLA-A*02:01)&P61769(B2M)&AAPGAPAAA
Molecular Weight:	35.6 kDa (HLA-A*02:01) and 11.9 kDa (B2M)

QC Testing

Biological Activity:	Immobilized Human HLA-A*02:01&B2M&PLK1 (AAPGAPAAA) Monomer, His Tag at 0.5 µg/ml (100 µl/well) on the plate. Dose response curve for Anti-B2M Antibody, mFc Tag with the EC50 of 12.7 ng/ml determined by ELISA.
Purity:	> 95% as determined by Bis-Tris PAGE; > 95% as determined by HPLC
Endotoxin:	< 1 EU/µg of the protein as determined by the LAL method.
Formulation:	Supplied as 0.22 µm filtered solution in 20 mM Tris, 200 mM NaCl (pH 8.0).

Preparation and Storage

Stability & Storage:

It is recommended to store recombinant proteins at -20°C to -80°C for future use. Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Protein Background

A key regulator of cell division is polo-like kinase 1 (PLK1), a member of the serine/threonine kinase family that plays an important role in regulating the mitotic and meiotic cell cycle. The phosphorylation of specific substrates mediated by PLK1 controls nuclear envelope breakdown (NEBD), centrosome maturation, proper spindle assembly, chromosome segregation, and cytokinesis.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481