

## HLA-A\*03:01&B2M&IGA2 (FTWTPSSGK) Monomer Protein, Human, MHC (His & Avi)

### General Information

Synonyms:	FTWTPSSGK;IGA2
Protein Construction:	Gly25-Thr305 (HLA-A*03:01), Ile21-Met119 (B2M) and FTWTPSSGK peptide
Species:	Human
Expression Host:	E. coli
Accession:	NP_002107.3 (HLA-A*03:01)&P61769 (B2M)&FTWTPSSGK
Molecular Weight:	35.26 kDa (HLA-A*03:01) and 11.9 kDa (B2M)

### QC Testing

Biological Activity:	Immobilized Human HLA*03:01&B2M&IGA2 (FTWTPSSGK) Monomer, His Tag at 1 µg/ml (100 µl/well) on the plate. Dose response curve for Anti-B2M Antibody, mFc Tag with the EC50 of 2.4 ng/ml determined by ELISA (QC Test).
Purity:	> 95% as determined by Bis-Tris PAGE; > 95% as determined by HPLC
Endotoxin:	< 1 EU/µg of the protein as determined by the LAL method.
Formulation:	Supplied as 0.22 µm filtered solution in 20 mM Tris, 200 mM NaCl (pH 8.0).

### Preparation and Storage

#### Stability & Storage:

It is recommended to store recombinant proteins at -20°C to -80°C for future use. Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

#### Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

### Protein Background

The immunoglobulin variable region is a part of an antibody molecule responsible for antigen recognition and binding. It is located at the N-terminal end of each antibody chain and is highly diverse due to the presence of variable (V), diversity (D), and joining (J) gene segments in B cells that undergo somatic recombination during antibody production. This variability allows the immune system to recognize a wide range of antigens.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

**This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use**

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481