

HLA-A*02:01&B2M&Foxp3 (TLIRWAILEA) Monomer Protein, Human, MHC (His & Avi), Biotinylated

General Information

Synonyms: TLIRWAILEA

Protein Construction: Gly25-Thr305 (HLA-A*02:01), Ile21-Met119 (B2M) and TLIRWAILEA peptide

Species: Human

Expression Host: E. coli

Accession: A0A140T913(HLA-A*02:01)&P61769(B2M)&TLIRWAILEA

Molecular Weight: 35.6 kDa (HLA-A*02:01) and 11.9 kDa (B2M)

QC Testing

Biological Activity: Immobilized Biotinylated Human HLA-A*02:01&B2M&Foxp3 (TLIRWAILEA) Monomer, His-Avi Tag at 1 µg/ml (100 µl/well) on the streptavidin precoated plate (5 µg/ml). Dose response curve for Anti-B2M Antibody, mFc Tag with the EC50 of 19.7 ng/ml determined by ELISA.

Purity: > 95% as determined by Bis-Tris PAGE; > 95% as determined by HPLC

Endotoxin: < 1 EU/µg of the protein as determined by the LAL method.

Formulation: Supplied as 0.22 µm filtered solution in 20 mM Tris, 200 mM NaCl (pH 8.0).

Preparation and Storage**Stability & Storage:**

It is recommended to store recombinant proteins at -20°C to -80°C for future use. Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Protein Background

TLIRWAILEA is position 344-353 of Foxp3 which is a member of a large forkhead family of related proteins. The peptide has minimal homology with other Foxp family members, such as Foxp1, 2 and 4, and it has also been shown to induce strong peptide-specific CD8+ T cell responses, which recognize Foxp3+/HLA-A*02:01+ cutaneous T lymphoma cells.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481