

HLA-A*02:01&B2M&PRAME (ALYVDSLFFL) Monomer Protein, Human, MHC (His & Avi), Biotinylated

General Information

Synonyms: OIP4;OIP-4;PRAME;MAPE

Protein Construction: Gly25-Thr305 (HLA-A*02:01), Ile21-Met119 (B2M) and ALYVDSLFFL peptide

Species: Human

Expression Host: HEK293 Cells

Accession: P79603(HLA-A*02:01)&P61769(B2M)&ALYVDSLFFL

Molecular Weight: 50.60 kDa (Predicted); 50-65 kDa (Reducing condition, due to glycosylation)

QC Testing

Biological Activity: Immobilized Biotinylated Human HLA-A*02:01&B2M&PRAME (ALYVDSLFFL) Monomer, His Tag at 1 µg/ml (100 µl/well) on the streptavidin precoated plate (5 µg/ml). Dose response curve for Anti-HLA-A*02:01&B2M&PRAME (ALYVDSLFFL) Antibody, hFc Tag with the EC50 of 6.2 ng/ml determined by ELISA.

Purity: > 95% as determined by Bis-Tris PAGE; > 95% as determined by HPLC

Endotoxin: < 1 EU/µg of the protein as determined by the LAL method.

Formulation: Lyophilized from 0.22 µm filtered solution in PBS, 100 mM L-arginine (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.

Preparation and Storage**Reconstitution:**

Reconstitute the lyophilized protein in distilled water. The product concentration should not be less than 100 µg/ml. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.

Stability & Storage:

It is recommended to store recombinant proteins at -20°C to -80°C for future use. Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Protein Background

PRAME (PReferentially expressed Antigen in MELanoma) is a melanoma-associated antigen expressed in cutaneous and ocular melanomas and some other malignant neoplasms, while its expression in normal tissue and benign tumors is limited.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481