

HLA-A*02:01&B2M&MAGE-A10 (GLYDGMEHL) Tetramer Protein, Human, MHC (His & Avi), PE-

General Information

Synonyms: MAGE-A10

Protein Construction: Gly25-Thr305 (HLA-A*02:01), Ile21-Met119 (B2M) and GLYDGMEHL peptide

Species: Human

Expression Host: HEK293 Cells

Accession: A0A140T913(HLA-A*02:01)&P61769(B2M)&GLYDGMEHL

QC Testing

Biological Activity: Activity has not been tested. It is theoretically active, but we cannot guarantee it.

Formulation: Lyophilized from 0.22 µm filtered solution in PBS, 0.2% BSA (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.

Preparation and Storage

Reconstitution:

Reconstitute the lyophilized protein in distilled water. The product concentration should not be less than 100 µg/ml. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.

Stability & Storage:

It is recommended to store recombinant proteins at -20°C to -80°C for future use. Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Protein Background

MAGE-A10 is a subtype of the Melanoma-associated antigen A (MAGE-A), a class of tumor antigens that are extensively expressed in various histological types of tumors and represents an attractive target for tumor immunotherapy. High-level expression of MAGE-A10 improved the anti-tumor immune cytotoxicity of MAGE-A10-specific CTLs in lung cancer cell lines and primary lung cancer cells.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481