

RNF43 Protein, Canine, Recombinant (His)

General Information

Synonyms:	RNF43;E3 ubiquitin-protein ligase RNF43;RING-type E3 ubiquitin transferase RNF43;RING finger protein 43
Protein Construction:	Gly24-Tyr197
Species:	Canine
Expression Host:	HEK293 Cells
Accession:	XP_038533015.1
Molecular Weight:	20.68 kDa (Predicted); 30-40 kDa (Reducing condition, due to glycosylation)

QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it.
Purity:	> 95% as determined by Bis-Tris PAGE; > 95% as determined by HPLC
Endotoxin:	< 1 EU/μg of the protein as determined by the LAL method.
Formulation:	Lyophilized from 0.22 μm filtered solution in PBS (pH 7.4); Normally 8% trehalose is added as protectant before lyophilization.

Preparation and Storage

Reconstitution:

Reconstitute the lyophilized protein in distilled water. The product concentration should not be less than 100 μg/ml. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.

Stability & Storage:

It is recommended to store recombinant proteins at -20°C to -80°C for future use. Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Protein Background

RNF43 (E3 ubiquitin-protein ligase RNF43 or RING-type E3 ubiquitin transferase RNF43) functions as a tumor suppressor, by exerting a predominant negative feedback mechanism in the Wnt/β-catenin signaling pathway. RNF43 inhibits Wnt/beta-catenin signaling by ubiquitinating Frizzled receptor and targeting it to the lysosomal pathway for degradation.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481