

SOST Protein, Mouse, Recombinant (His & Avi), Biotinylated

General Information

Protein Construction: Gln24-Tyr211

Species: Mouse

Expression Host: HEK293 Cells

Accession: Q99P68

Molecular Weight: 24.04 kDa (Predicted); 35-50 kDa (Reducing condition, due to glycosylation)

QC Testing

Biological Activity: Immobilized Biotinylated Mouse SOST, His Avi Tag at 1 µg/ml (100 µl/well) on the streptavidin precoated plate (5 µg/ml). Dose response curve for Human LRP-6, mFc Tag with the EC50 of 13.2 ng/ml determined by ELISA (QC Test).

Purity: > 95% as determined by Bis-Tris PAGE

Endotoxin: < 1 EU/µg of the protein as determined by the LAL method.

Formulation: Supplied as 0.22 µm filtered solution in PBS (pH 7.4).

Preparation and Storage

Stability & Storage:

It is recommended to store recombinant proteins at -20°C to -80°C for future use. Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Protein Background

SOST, also known as sclerostin, is a member of the cerberus/DAN family, a group of secreted glycoproteins characterized by a cysteine-knot motif. SOST is negative regulator of bone growth that acts through inhibition of Wnt signaling and bone formation.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481