

HLA-A*03:01&B2M&KRAS G12V (VVGAVGVGK) Monomer Protein, Human, MHC (His & Avi),

General Information

Synonyms: KRAS2;K-RAS4A;GTPase Kras;KI-RAS;CFC2;K-Ras 2;NS;MHC;RALD;C-K-RAS;RASK2;KRAS1;K-RAS4B;NS3;KRAS;K-RAS2A;K-RAS2B

Protein Construction: Gly25-Thr305 (HLA-A*03:01), Ile21-Met119 (B2M) and VVGAVGVGK peptide

Species: Human

Expression Host: HEK293 Cells

Accession: NP_002107.3(HLA-A*03:01)&P61769(B2M)&VVGAVGVGK

Molecular Weight: 50.09 kDa (Predicted); 51-60 kDa (Reducing condition, due to glycosylation)

QC Testing

Biological Activity: Activity has not been tested. It is theoretically active, but we cannot guarantee it.

Purity: > 95% as determined by Bis-Tris PAGE; > 95% as determined by HPLC

Endotoxin: < 1 EU/μg of the protein as determined by the LAL method.

Formulation: Lyophilized from 0.22 μm filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.

Preparation and Storage

Reconstitution:

Reconstitute the lyophilized protein in distilled water. The product concentration should not be less than 100 μg/ml. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.

Stability & Storage:

It is recommended to store recombinant proteins at -20°C to -80°C for future use. Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Protein Background

Kirsten rat sarcoma 2 viral oncogene homolog (KRAS) is the most commonly mutated oncogene in human cancer. The developments of many cancers depend on sustained expression and signaling of KRAS, which makes KRAS a high-priority therapeutic target. The virtual screening approach to discover novel KRAS inhibitors and synthetic lethality interactors of KRAS are discussed in detail.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481