

2'-OMe-U(Bz)-CE-Phosphoramidite

Chemical Properties

CAS No. :	178422-03-2
Formula:	C47H53N4O10P
Molecular Weight:	864.93
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	2'-OMe-U(Bz)-CE-Phosphoramidite is K00145-CP, used for oligonucleotide synthesis research.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.1562 mL	5.7808 mL	11.5616 mL
5 mM	0.2312 mL	1.1562 mL	2.3123 mL
10 mM	0.1156 mL	0.5781 mL	1.1562 mL
50 mM	0.0231 mL	0.1156 mL	0.2312 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481