

5'(E)-VP-2'-OMe-U-CE-Phosphoramidite

Chemical Properties

CAS No. :

Formula: C₃₂H₅₂N₄O₁₃P₂

Molecular Weight: 762.73

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	5'(E)-VP-2'-OMe-U-CE-Phosphoramidite is VP OMe-U monomer, used for oligonucleotide synthesis research.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.3111 mL	6.5554 mL	13.1108 mL
5 mM	0.2622 mL	1.3111 mL	2.6222 mL
10 mM	0.1311 mL	0.6555 mL	1.3111 mL
50 mM	0.0262 mL	0.1311 mL	0.2622 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481