

Alkyne-Modifier SerinolPhosphoramidite

Chemical Properties

CAS No. :	2084118-98-7
Formula:	C42H55N4O7P
Molecular Weight:	758.9
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	Alkyne-Modifier SerinolPhosphoramidite is Spacer monomer, used for oligonucleotide synthesis research.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.3177 mL	6.5885 mL	13.177 mL
5 mM	0.2635 mL	1.3177 mL	2.6354 mL
10 mM	0.1318 mL	0.6588 mL	1.3177 mL
50 mM	0.0264 mL	0.1318 mL	0.2635 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481