

2'-TC-U -3'-CE-Phosphoramidite

Chemical Properties

CAS No. :	1219089-92-5
Formula:	C44H54N5O11PS2
Molecular Weight:	924.03
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	2'-TC-U -3'-CE-Phosphoramidite can be used in research related to oligonucleotide synthesis.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.0822 mL	5.4111 mL	10.8222 mL
5 mM	0.2164 mL	1.0822 mL	2.1644 mL
10 mM	0.1082 mL	0.5411 mL	1.0822 mL
50 mM	0.0216 mL	0.1082 mL	0.2164 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481