

THA8 acid

Chemical Properties

CAS No. :

Formula: C78H125N7O36

Molecular Weight: 1736.87

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	THA8 acid can be used in research related to oligonucleotide synthesis.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.5757 mL	2.8787 mL	5.7575 mL
5 mM	0.1151 mL	0.5757 mL	1.1515 mL
10 mM	0.0576 mL	0.2879 mL	0.5757 mL
50 mM	0.0115 mL	0.0576 mL	0.1151 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481