

## TA14 Amidite

## Chemical Properties

CAS No. :

Formula: C51H86N7O13P

Molecular Weight: 1036.26

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	TA14 Amidite can be used in research related to oligonucleotide synthesis.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.965 mL	4.825 mL	9.6501 mL
5 mM	0.193 mL	0.965 mL	1.930 mL
10 mM	0.0965 mL	0.4825 mL	0.965 mL
50 mM	0.0193 mL	0.0965 mL	0.193 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481