

## Bis-(m-PEG4)-amidohexanoic acid

## Chemical Properties

CAS No. : 2353409-74-0


Formula: C<sub>26</sub>H<sub>50</sub>N<sub>2</sub>O<sub>12</sub>

Molecular Weight: 582.68

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Bis-(m-PEG4)-amidohexanoic acid is a polyethylene glycol (PEG)-based linker utilized in the synthesis of proteolysis-targeting chimeras (PROTACs)[1].
Targets(IC50)	Others,PROTAC Linker
In vitro	PROTACs consist of two ligands linked together: one binds to an E3 ubiquitin ligase and the other to the target protein. These compounds utilize the intracellular ubiquitin-proteasome system to selectively degrade target proteins[1].

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.7162 mL	8.581 mL	17.1621 mL
5 mM	0.3432 mL	1.7162 mL	3.4324 mL
10 mM	0.1716 mL	0.8581 mL	1.7162 mL
50 mM	0.0343 mL	0.1716 mL	0.3432 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

## Reference

An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481