


G-749

Chemical Properties

CAS No. :	1457983-28-6
Formula:	C ₂₅ H ₂₅ BrN ₆ O ₂
Molecular Weight:	521.41
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	G-749, a new-type FLT3 inhibitor, exhibits effective and sustained inhibition of the FLT3 wild type and mutants.
Targets(IC50)	Apoptosis,FLT,c-RET,Aurora Kinase,TAM Receptor
In vitro	In primary cells from AML patients, G-749 was able to inhibit cell proliferation by inducing apoptosis. In BaF3 cell lines stably expressing FLT3-ITD/N676D, FLT3-ITD/F691L, FLT3-D835Y, or FLT3-D835Y/N676D (IC ₅₀ <10 nM), G-749 showed significant inhibitory activity.
In vivo	In primary cells from AML patients, G-749 was able to inhibit cell proliferation by inducing apoptosis. In BaF3 cell lines stably expressing FLT3-ITD/N676D, FLT3-ITD/F691L, FLT3-D835Y, or FLT3-D835Y/N676D (IC ₅₀ <10 nM), G-749 showed significant inhibitory activity.
Kinase Assay	Kinase assay: Activity assays are conducted using Lance Ultra time-resolved fluorescence resonance energy transfer (TR-FRET) technology from Perkin-Elmer. Briefly, 10 ng/mL FLT3 enzyme, a serial diluted G-749, 80 nM substrate of ULight-poly-GT peptide and variable amounts of ATP (8.5 μM to 1088 μM) are mixed in kinase assay buffer (50 mM HEPES pH 7.5, 10 mM MgCl ₂ , 1 mM EGTA, 2 mM DTT and 0.01% Tween-20) and are added to a 384-well OptiPlate-384 in a volume of 10 μL. Kinase reactions are incubated at room temperature for up to 1 h and then stopped by the addition of 5 μL of 10 mM EDTA. A volume of 5 μL of the specific Eu-labeled-anti-phosphopeptide antibody diluted in LANCE Detection Buffer is then added to a final concentration of 2 nM. After 30-minute incubation, assay plates are incubated at 23°C and the LANCE signal is measured on an EnVision Multilabel Reader. Excitation wavelength is set at 320 nm and emission monitored at 615 nm (donor) and 665 nm (acceptor). The IC ₅₀ is calculated using nonlinear regression analysis analysis by GradPad Prism 5.
Cell Research	Cells are seeded at a density of 2 ×10 ⁴ cells per well and treated with the indicated concentrations of test inhibitor for 72 hours at 37°C. The conditioned medium (CM) is prepared from HS-5 cell culture for 5 days under routine culture conditions, clarified by centrifugation, and used immediately. The CM is added to complete medium at a final concentration of 35%. In coculture experiments, 5 ×10 ⁴ AML blast cells are plated in 24-well plates containing 1 ×10 ⁴ HS-5 monolayers and then cultured for at least 48 hours before the exposure of inhibitors. Cell viability is determined by an ATPLite assay. (Only for Reference)

Solubility Information

Solubility	Ethanol: < 1 mg/mL (insoluble or slightly soluble), H2O: < 1 mg/mL (insoluble or slightly soluble), DMSO: < 1 mg/mL (insoluble or slightly soluble), (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9179 mL	9.5894 mL	19.1788 mL
5 mM	0.3836 mL	1.9179 mL	3.8358 mL
10 mM	0.1918 mL	0.9589 mL	1.9179 mL
50 mM	0.0384 mL	0.1918 mL	0.3836 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Lee HK, et al. Blood. 2014, 123(14), 2209-2219.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481