


Phenol, m-(1-(p-hydroxyphenethyl)-3-isopropyl-3-pyrrolidinyl)-, hydrochloride, hydrate

Chemical Properties

CAS No. :	64047-91-2
Formula:	C ₂₁ H ₃₀ ClNO ₃
Molecular Weight:	379.92
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Phenol, m-(1-(p-hydroxyphenethyl)-3-isopropyl-3-pyrrolidinyl)-, hydrochloride, hydrate is a bioactive chemical.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.6321 mL	13.1607 mL	26.3213 mL
5 mM	0.5264 mL	2.6321 mL	5.2643 mL
10 mM	0.2632 mL	1.3161 mL	2.6321 mL
50 mM	0.0526 mL	0.2632 mL	0.5264 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481