


WEN05-03

Chemical Properties

CAS No. : 2419180-34-8
 Formula: C₁₇H₁₄N₄O₄S
 Molecular Weight: 370.38
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	WEN05-03 is a novel inhibitor of type III secretion system ATPase EscN from enteropathogenic E. coli.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.6999 mL	13.4996 mL	26.9993 mL
5 mM	0.540 mL	2.6999 mL	5.3999 mL
10 mM	0.270 mL	1.350 mL	2.6999 mL
50 mM	0.054 mL	0.270 mL	0.540 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481