


Stachartin C

Chemical Properties

CAS No. : 1978388-56-5
 Formula: C₂₉H₄₁NO₆
 Molecular Weight: 499.64
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	Stachartin C is a natural product from the cultures of the tin mine tailings-associated fungus <i>Stachybotrys chartarum</i> .
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0014 mL	10.0072 mL	20.0144 mL
5 mM	0.4003 mL	2.0014 mL	4.0029 mL
10 mM	0.2001 mL	1.0007 mL	2.0014 mL
50 mM	0.040 mL	0.2001 mL	0.4003 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481