

## 4-Phenoxy pyridine

## Chemical Properties

CAS No. :	4783-86-2
Formula:	C <sub>11</sub> H <sub>9</sub> NO
Molecular Weight:	171.195
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

## Biological Description

Description	4-Phenoxy pyridine ,with CAS No. 4783-86-2, is a fragment molecule that serves as an important scaffold for molecular linking, expansion, and modification. 4-Phenoxy pyridine provides a structural basis and research tool for the design and screening of novel drug candidates, and is commonly used in drug discovery, drug synthesis, and related research.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.8411 mL	29.2056 mL	58.4112 mL
5 mM	1.1682 mL	5.8411 mL	11.6822 mL
10 mM	0.5841 mL	2.9206 mL	5.8411 mL
50 mM	0.1168 mL	0.5841 mL	1.1682 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481