


Linderane

Chemical Properties

CAS No. :	13476-25-0
Formula:	C ₁₅ H ₁₆ O ₄
Molecular Weight:	260.29
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	1. Linderane was characterized as a mechanism-based inactivator of CYP2C9.
Targets(IC50)	Cytochromes P450

Solubility Information

Solubility	DMSO: 47 mg/mL (180.57 mM), Sonication is recommended. Chloroform, Dichloromethane, Ethyl Acetate, Acetone, etc.: Soluble, (< 1 mg/ml refers to the product slightly soluble or insoluble)
In vivo Formulation	10% DMSO+90% Corn Oil: 2 mg/mL (7.68 mM), Sonication is recommended. <i>Please add the solvents sequentially, clarifying the solution as much as possible before adding the next one. Dissolve by heating and/or sonication if necessary. Working solution is recommended to be prepared and used immediately. The formulation provided above is for reference purposes only. In vivo formulations may vary and should be modified based on specific experimental conditions.</i>

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.8419 mL	19.2093 mL	38.4187 mL
5 mM	0.7684 mL	3.8419 mL	7.6837 mL
10 mM	0.3842 mL	1.9209 mL	3.8419 mL
50 mM	0.0768 mL	0.3842 mL	0.7684 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Determination of linderane in root tuber of *Lindera aggregata* by HPLC,

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481