

Ald-Ph-PEG4-NH-Boc

Chemical Properties

CAS No. : 2110449-01-7


Formula: C₂₁H₃₂N₂O₈

Molecular Weight: 440.49

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Ald-Ph-PEG4-NH-Boc, a PEG-based linker for PROTACs, connects two essential ligands, crucial for forming PROTAC molecules. This linker facilitates selective protein degradation by leveraging the ubiquitin-proteasome system within cells.
Targets(IC50)	Others,PROTAC Linker
In vitro	PROTACs consist of two distinct ligands connected by a linker: one ligand targets an E3 ubiquitin ligase, and the other targets the desired protein. They utilize the intracellular ubiquitin-proteasome system to selectively degrade these target proteins[1].

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2702 mL	11.351 mL	22.702 mL
5 mM	0.454 mL	2.2702 mL	4.5404 mL
10 mM	0.227 mL	1.1351 mL	2.2702 mL
50 mM	0.0454 mL	0.227 mL	0.454 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481